



MennMove 300
Three-channel Electrocardiograph



www.mennenmedical.com

Israel: Mennen Medical Ltd. 4 Ha-Yarden Street, Yavne, P.O. Box 102, Rehovot 7610002, Israel. Phone: +972-8-9323333. Fax: +972-8-9328510
United States: Mennen Medical Corp. 290 Andrews Road, Feasterville-Trevose, PA 19053-3480, USA. Phone: +215 259-1020. Fax: +215 357-2010
United Kingdom: Charter Kontron Ltd. Unit 18 Avant Business Center 21 Denbigh Road Milton Keynes MK1 1DT, UK. Phone: +44 1908 646070. Fax: +44 1908 646030



We Follow Every Beat of Your Heart

MennMove 300

Three-channel Electrocardiograph



Mennen Medical Ltd. All rights reserved. Data is subject to be changed without notice.

MennMove 300

Three-channel Electrocardiograph



⏪ MennMove 300

MennMove 300 Features:

Portable & User-friendly Design

- Vivid waveforms display
- Built-in rechargeable Li-ion battery, AC/DC power supply
- Communication interface: Ethernet/RS232/USB

Streamlined Workflow

- One-button operation
- Auto/Manual/Rhythm work modes
- Automatic baseline adjustment
- Lead-off detection and alarm

Accurate Clinical Performances

- 12 leads simultaneous acquisition with optional 24-bit A/D converter
- Automatic ECG measurements and interpretation tested with CSE/AHA/MIT database
- Complete digital filtering, resisting baseline drift, AC and EMG interference

Flexible Printing Solutions


- Built-in high-resolution thermal printer
- 12 leads grid printing on normal A4 paper by external printer

Data Management

- Large patient data storage: 144 ECGs, which can be enlarged by a USB flash disk
- PC-based ECG data management software (optional)



PC-based ECG data management software



Smart ECG Viewer

Intuitive interface to improve your work efficiency

