

Benefits

- Bidirectional communication with Hospital Information Systems in various protocols: HL7, DICOM, XML etc*.
- Interface to Mennen's XVu systems - Automatic inclusion of patient's monitored data and notes made during procedure in patient record.
- Physician report builder for unlimited customization of reports.
- Unlimited database query capability
- Hemodynamic calculations and waveforms
- Locally acquired images available for inclusion in the patient record
- Coronary artery diagnostic and intervention findings
- Designated field for free text summary and recommendations
- Incorporation of all Peripheral Diagrams
- Flexible statistical report generator
- Inventory and Billing package
- Scheduling module
- Interface to Mennen's Patient Monitors (VitaLogik & Envoy)*

*Available with server only

Pediatric Package - includes all pediatric cardiologist necessities

- Mullins pediatric diagrams
- Complete pediatric calculations
- Easy access to historical data
- Drag and drop of Saturations and Values

Mennen Medical has wide experience and knowledge with customizing systems connectivity, providing each hospital with its specific integration requirements and unique properties.



www.mennenmedical.com



We Follow Every Beat of Your Heart

HemoCIS plus

Cardiology Information System



Headquarters:
Mennen Medical Ltd.
4 Ha-Yarden Street,
Yavne, P.O. Box 102,
Rehovot 76100, Israel
Phone: +972-8-9323333
Fax: +972-8-9328510



United States:
Mennen Medical Corp.
290 Andrews Road,
Feasterville-Treose,
PA 19053-3480, USA
Phone: +215 259-1020
Fax: +215 357-2010

EC REP



United Kingdom:
Charter Kontron Ltd.
Unit 18 Avant Business Center
21 Denbigh Road, Milton
Keynes MK1 1DT, UK
Phone: +44 1908 646070
Fax: +44 1908 646030

Enhancing Data Integration

The **HemoCIS plus** is a powerful database that is capable of storing a large volume of patient data. In addition, it provides easy, rapid and efficient data access to patient medical report system.

The **HemoCIS plus** includes in a single system, many necessities of Cardiac Catheterization wards. For example:

- Scheduling
- Inventory
- Billing
- Patient files
- Outcome Reporting
- Statistics




Assigning user rights:

Each user is authorized to either view or modify specific fields.

The **HemoCIS plus** with a server provides several bidirectional CONNECTIVITY capabilities. This system enables users to connect several Horizon XVu cardiac procedure rooms to:

- HIS ADT by HL7
- DICOM worklist
- Laboratory Information System (LIS)
- Interface to PACS
- Hospital Billing System
- Reporting by HL7 and XML


Modern Graphic User Interface (GUI):

<p>Unlimited customization of reports</p>	<p>Patient file management</p>	<p>"Outlook" like Scheduler</p>
		


Connect to Horizon XVu

The HemoCIS plus can be installed in two optional configurations:

1. **HemoCIS plus on the Horizon XVu**- designed with the ability to run on the Horizon XVu PC in a cardiac unit with one procedure room as an integrated system. The staff can utilize the features of the hemodynamic system and HemoCIS plus on one combined station.



2. **HemoCIS plus with a Server** - Available to users with multiple procedure rooms, enables the use of multiple HemoCIS clients connected via a HemoCIS server.



From admission until discharge

<p>Patient Registration</p> <ul style="list-style-type: none"> • Demographic Data • Procedure Scheduling • DICOM Worklist* 	<p>Review</p> <ul style="list-style-type: none"> • XVu Hemodynamic data and Procedure Notes 	<p>Findings</p> <ul style="list-style-type: none"> • Inventory • Statistical Analysis • Coronary arteries diagnostic and intervention 	<p>Final Report</p> <ul style="list-style-type: none"> • Generate Report • Send final reports to: <ul style="list-style-type: none"> • HIS* • Billing* • Hospital Reporting system* • Hard copy print
--	---	---	---

*Option