

EP Combo System

Horizon XVu + EMS XL Cardiovascular Workstation

The Mennen Medical Combo system combines on the same consol the EMS-XL cardiac electrophysiology and the HORIZON XVu Computerized Cathlab Hemodynamic System.

The Combo provides a smooth solution for hospitals that use the same procedure room for Cardiac Electrophysiology and Hemodynamic catheterization.

Features

- A Video Switch enables the user to switch the system from Electrophysiology mode to Hemodynamic mode without need to replace or change display screens, keyboard or mouse.
- The same network printer serves both applications.
- The front end amplifier for the Horizon XVu is the CFE.
- The front end amplifier for the EMS-XL is a 32 or 64 channel surface and intra-cardiac ECG amplifier with built in dual channel stimulator.

EMS-XL



Horizon XVu



For more information:

Horizon XVu Brochure for Hemodynamic system details

EMS-XL Brochure for Cardiac electrophysiology detail

www.mennenmedical.com