

CritiCool MINI

COMPACT design, COMPLETE treatment



Battery operated
thermoregulation device

Fits ground & aerial
transportation



Unique Core
Temperature Splitter

Advanced & intuitive
user interface



A new **SMALL & PORTABLE** device

Designed for a **QUICK & EASY** set up

- Color Touch screen
- Designed for ground and aerial transport
- Fits different shelving options in different hospital department
- One hour battery operation span
- Built-in unique Core Temperature Splitter
- 10 day data storage
- Built-in data viewer
- A new Thermal Disinfection application

Fast

- Optimal solution for early initiation of TTM

Simple

- Friendly graphic user interface
- Small system, easy transport

Precise Control

- Precise TTM, programmed Set Point temperature
- Automatic and stable temperature regulation and monitoring

Set

- Intuitive and simple settings of temperature

Wrap

- Apply the CureWrap garment

Cool & Rewarm

- Cool fast
- Step wise rewarming



* Cart optional

Israel: MTRE Advanced Technologies Ltd.
6 Ha-Kishon Street, Yavne.
P.O. Box 102, Rehovot 7610002, Israel.
Phone: +972-8-9323333.
Fax: +972-8-9328510
www.mtre.com

United States: Mennen Medical Corp.
290 Andrews Road, Feasterville-Trevoze, PA
19053-3480, USA.
Phone: +215 259-1020.
Fax: +215 357-2010
www.mennenmedical.com

United Kingdom: Charter Kontron Ltd.
Unit 18 Avant Business Center 21
Denbigh Road Milton Keynes MK1 1DT, UK.
Phone: +44 1908 646070.
Fax: +44 1908 646030
www.charter-kontron.com